Genarrayt® Microarray

Genarrayt™

Genesis Diagnostics Food Intolerance Microarray Technology

Beta Diagnostici sas
Viale Europa is 57 n°22 98123 MESSINA
Tel. 090 696001 fax 090 695563
biotechnology@betadiagnostici.it - info@betadiagnostici.it
www.betadiagnostici.it


 **BETA DIAGNOSTICI**
biotechnology

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

Methods

Results

Proprietary software converts raw data files into individual patient reports


Results are then discussed between the patient and a nutritionist who will suggest dietary changes if applicable

FOODPRINT 2000+

Completed report for:
Patient Number: 0
Sample Date: 03-Mar-07

Food	Result	Food	Result	Food	Result	Food	Result
APPLE	negative	APPLE	negative	APPLE	negative	APPLE	negative
BANANA	negative	BANANA	negative	BANANA	negative	BANANA	negative
BREAD	negative	BREAD	negative	BREAD	negative	BREAD	negative
BUTTER	negative	BUTTER	negative	BUTTER	negative	BUTTER	negative
CHEESE	negative	CHEESE	negative	CHEESE	negative	CHEESE	negative
EGG	negative	EGG	negative	EGG	negative	EGG	negative
FISH	negative	FISH	negative	FISH	negative	FISH	negative
GRAPE	negative	GRAPE	negative	GRAPE	negative	GRAPE	negative
ORANGE	negative	ORANGE	negative	ORANGE	negative	ORANGE	negative
PANCAKE	negative	PANCAKE	negative	PANCAKE	negative	PANCAKE	negative
PASTA	negative	PASTA	negative	PASTA	negative	PASTA	negative
PEPPER	negative	PEPPER	negative	PEPPER	negative	PEPPER	negative
PIZZA	negative	PIZZA	negative	PIZZA	negative	PIZZA	negative
RICE	negative	RICE	negative	RICE	negative	RICE	negative
SOY	negative	SOY	negative	SOY	negative	SOY	negative
TOMATO	negative	TOMATO	negative	TOMATO	negative	TOMATO	negative
WATERMELON	negative	WATERMELON	negative	WATERMELON	negative	WATERMELON	negative
YOGURT	negative	YOGURT	negative	YOGURT	negative	YOGURT	negative

BRINGING INNOVATION TO IMMUNOASSAY

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




Options...

Help

Image
 Preview (fast scan)
 Wavelength 532

Tools


79 73


Feature Viewer
Block: Feature:
X, Y (µm):
Name:

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
Cutoff at IgG concentration (U/mL):

Concentration at borderline (U/mL):

8759	3940	8821	4702	5932	Delete outliers		
5329	5626	5880	7521		Reset with 'New Patient' button		

Mean control values	
Positive	Negative
94	4

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Patient No: 995678 Test Ref: 678

Test date: 10-Jul-08

Buttons: Prepare ScanData (Prepared), New Patient, Print This Pad Data, Save Excel File, Compile Data (Data Compiled)

17	negative	BEAN (white ha	negative	LEEK	negative	OX	negative	CHERRY	negative	POMEGRANAT	negative	CUTTLEFISH	negative	PLAICE
16	negative	BEAN (red kidn	negative	LENTIL	negative	HORSE	negative	COCONUT	negative	RASPBERRY	negative	DORADO	negative	RAZOR CLAM
15	negative	BEETROOT	negative	LETT	negative	BILLY GOAT	negative	DATE	negative	REDCURRANT	negative	EEL	negative	SALMON
14	negative	BROCCOLI	negative	MAR	negative	PORK	negative	FIG	negative	STRAWBERRY	POS	POS	POS	
13	negative	CABBAGE	negative	MUS	negative	QUAIL	negative	GRAPE	negative	WATERMELON	0	0	0	
12	negative	CARROT	negative	ONIO	negative	RABBIT	negative	GRAPEFRU	negative	AGAR AGAR	40	40	40	
11	negative	CAULIFLOWER	negative	PAT	negative	LAMB	negative	GUAYA	negative	ALGA espaguet	160	160	160	
10	negative	CAYENNE	negative	PEA	negative	WILD BOAR	negative	KIWI	negative	ALGA spirulina				
9	negative	CELERY	negative	POTATO	negative	DUCK	negative	LEMON	negative	ALGA wakame	negative	CRANBERRY		
8	negative	CHARD	negative	GOURD (squas	negative	TURKEY	negative	LIME	negative	ANCHOVY	negative	HOPS		
7	negative	CHICKPEA	negative	RADISH	negative	PARTRIDGE	negative	LYCHEE	negative	BARNACLE	negative	HADDOCK		
6	negative	CHICORY	negative	CABBAGE (red	negative	CHICKEN	negative	MANGO	negative	BASS	negative	SCALLOP		
5	negative	CUCUMBER	negative	ROCKET	negative	YEAL	negative	MELON (Ho	negative	CARP	negative	SHAKE	negative	SARDINE
4	negative	DILL	negative	SPINACH	negative	BEEF	negative	MULBERRY	negative	AVOCADO	negative	HERRING	negative	SHRIMP/PRA
3	negative	FENNEL (leaf)	negative	BRUSSEL SPF	negative	VENISON	negative	NECTARINE	negative	TANGERINE	negative	LOBSTER	negative	SOLE
2	51	Negative contr	negative	SWEET POTAT	negative	APPLE	negative	OLIVE	negative	WALNUT	negative	MACKEREL	negative	SQUID
1	173	Positive contr	negative	TAPIOCA	negative	APRICOT	negative	ORANGE	negative	ASPARAGUS	negative	MONKFISH	negative	SWORDFISH

IgG Standard curve

Calibrator	OD 1	OD 2	OD 3	Mean OD	SD	CV%	OD - Zero
0	100	100	100	100			0
160		37370		37370	4523	12	37270


OD/Cutoff at IgG concentration: 30 = 6988



OD/concentration at borderline: 20 = 4659

47678	46644	45124	47258	42678	Delete outliers
38531	36884	37595	39162	47991	Reset with 'New Patient' button

Contol values

Positive	Negative
178	103
173	51
170	
137	
140	

FOODPRINT™ 40+

Compiled report for :

Patient Number: 0

Test Date: 00-Jan-00

0 40+

40+	Type	Avoid (>30)	Borderline (20-30)	No reaction
15 COW'S MILK	dairy			
16 EGG WHITE	dairy			
17 EGG YOLK	dairy			
31 COD	fish/seafood			
32 CRAB	fish/seafood			
36 HADDOCK	fish/seafood			
39 LOBSTER	fish/seafood			
47 PLAICE	fish/seafood			
52 SHRIMP/PRAWN	fish/seafood			
61 APPLE	fruit			
65 BLACKBERRY	fruit			
73 GRAPEFRUIT	fruit			
76 LEMON	fruit			
84 ORANGE	fruit			
87 PEAR	fruit			
92 RASPBERRY	fruit			
95 STRAWBERRY	fruit			
100 BARLEY	grain			
102 CORN (maize)	grain			

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