**EVOLUTION
OF THE
REVOLUTION**

DRG:HYBRiD•XL[®]

Fully automated
Random Access
Analyzer for
Immunoassays and
Clinical Chemistry

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DRG

DRG:HYBRID•XL®

Fully automated Random Access Analyzer for the simultaneous measurement of Immunoassays and Clinical Chemistry

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Innovation

DRG introduces the totally newly designed **DRG:HYBRID•XL** Analyzer.

This innovative and unique technology allows, for the first time, the simultaneous measurement of immunoassays and clinical chemistry parameters including turbidimetric tests in one sample.

The comprehensive assay menu will include routine and niche parameters.

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User Orientation

The integrated touchscreen monitor and a well-defined software with simple user navigation guarantees easy system operation.

Ready-to-use reagents are provided in special reagent cartridges.

Bi-directional software-connection to laboratory systems via ASTM Interface is possible.

Samples and reagent cartridges as well as the master curve are easily identified by barcodes. Required periodic recalibration is done by 2-point calibration.

DRG:HYBRID-XL offers different working modes: Random Access and STAT possible.

Sample pre-dilution (up to 1:1000) can be performed automatically in the analyzer in separate dilution cartridges. Reagents are stable for 12 months from production date.

DRG offers packages of 20, 40, and 100 tests for flexibility and customization for unique test requirements. Small package sizes avoid needless and expensive waste due to expired reagents.

The sample needle is designed for clot detection and liquid level detection.

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
Efficiency

DRG:HYBRID-XL is designed as bench-top unit with small footprint. It allows you to process your samples totally automatically and efficiently; even in small lab rooms.

The device is incredibly quiet and is developed to be very user and service friendly. The rinse and washing solutions can be changed with just one hand. The instrument allows immediate determination „measurement on the fly“, generating the first results in only 10 minutes (Clinical Chemistry), respectively 60 minutes (Immunoassay), depending on the test with very short hands on times.

The device is able to process primary tubes from 13 to 16 mm diameter or defined secondary tubes.

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DRG:HYBRID-XL[®]

Fully automated Random Access Analyzer
for Immunoassays and **Clinical Chemistry**

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Specifications

Test Principle:

Enzyme Immunoassays
Clinical Chemistry
Immunoturbidimetry

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Samples/Cartridges:

Up to 20 samples/run
(12 primary tubes or 20 secondary tubes)
Random access possible
STAT possible
Reagents on board (cartridge)
(ready to use)
Barcoded reading of samples
and reagents
Serum, Plasma, Urine or Saliva samples
Sample dilution up to 1:1000

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Calibration:

Barcoded Master Curve
2 point recalibration
Long Stability

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Efficiency:

First results after 10 minutes (Clinical Chemistry)
and after 60 minutes (Immunoassay) respectively
True walk away capabilities
Short hands-on times
Benchtop Analyzer (60.8 x 63.5 x 58.6 cm high)
Only 2 system-reagents (Wash/Needle cleaning)
Clot detection
Liquid level detection
Variable assay-sizes 40 (20, 100 tests)

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Handling/Interface:

Color Monitor, Touchscreen (8.5 inch.)
Smart user guidance
Bi-directional Interface (ASTM)

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Flexibility

The **DRG:HYBRID•XL** fully automated analyzer is the ideal device for small and medium-sized labs. However, it is also suitable for high workload laboratories performing special routine and niche parameters totally automatically. Up to 40 different tests of up to 20 patient samples can be run. The analyzer is flexible enough to measure 40 different tests from one sample or just one or two tests from 20 samples.

Offering a wide range of adapted tests, this instrument can operate as a back-up instrument in case of breakdown of your large lab analyzer.

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Combination

The unique construction of the **DRG:HYBRID•XL** allows each user to combine individual diagnostic panels from immunoassays and clinical chemistry, e.g. a cardiac panel, a diabetes panel, or a metabolic syndrome panel.

DRG:HYBRID•XL can determine immunoassays quantitatively, as well as qualitatively (cut-off).

The analyzer is validated for detection of antigens and antibodies in serum, plasma or saliva.

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DRG Instruments GmbH, Germany
DRG International Inc. USA